



Lodox X-ray Technology: An Explanation¹

Conventional, Computed and Digital X-ray Systems




- Conventional, CR and flat panel detector-based DR systems use a wide cone-beam of X-rays which spread over the desired field-of-view in all directions (i.e. length and width).
- This wide beam results in a large number of scattered X-ray photons, which reduce the number of primary X-rays that reach the detector, and degrade the image quality by causing noise in the image.
- Therefore most systems employ a post-patient scatterrejection grid to reduce the number of scattered X-rays which would otherwise reach the detector or film.
- The result is an increase in overall patient radiation exposure in order to counteract the negative effects of scattered radiation from the cone-beam, and the absorption of radiation by the grid.
- This higher radiation exposure and the wide cone of X-rays also causes relatively high amounts of scattered radiation throughout the room.
- The geometric design of conventional, CR and DR systems limits the maximum possible field-of-view to approximately 400mm/16" square.

The Lodox Linear Slot Scanning Radiography System



- Lodox LSSR scanning system is a unique technology employing our proprietary X-ray beam controlling mechanism, and linear scanning technique.
- The Lodox X-ray beam is highly collimated by a narrow slit and a fan-width adjustment collimator into a laser-like fan beam of primary X-ray photons which spreads out in only one direction (i.e. width).
- The narrow beam and high quality primary photons reduce the number of X-rays scattered by the body, so no postpatient scatter-rejection grid is required.
- In addition, the detector is fitted into a scatter-absorbing housing, which is designed to eliminate almost all remaining scatter before detection, increasing the signal to noise ratio in the detected image.
- The inherently higher quality image means that a significantly lower patient exposure is required to achieve diagnostic-quality images.
- The lower exposure and collimated, narrow fan-beam also significantly reduce scattered radiation throughout the room.
- The X-ray tube, X-ray fan beam, collimating slit and detector all move together along a linear scanning path, collecting X-ray information to produce X-ray images of 100mm/4" square up to 1800mm/70" x 650mm/25".

Full-body, high-speed digital radiology with low radiation emission and scatter



Xmplar-dr Product Specifications

Image Quality >16000 grey levels; 60µm fundamental pixel size; up to 5 lp/mm

Image Size 1800mm x 680mm (70" x 26")

Throughput 35mm/s - 140mm/s; 12.98s for a full scan at normal speed

Duration of Examination <15 seconds from 'end of scan' until image available on screen

28 seconds between successive scans (provided X-ray tube heat <20%)

Positioning C-arm allows imaging angles from 0° - 90° i.e. AP (or PA) to lateral

Table can be tilted from 0 to 340mm, allowing Trendelenburg angles of 0° to 10°

Output Average instantaneous X-ray exposure time of 22ms

Maximum direct dose absorption is 1 mGy

Maximum scattered radiation measured 1 metre from the focal spot in any direction is 20uGy per hour

@ 130kV 25mA

Detector Proprietary ultra-low noise TDI CCD detector using Rarex Green Fast light conversion (Gadox)

Trolley/Gurney 220kg/440lbs weight capacity (with optional extra of 330kgs/660lbs bariatric trolley)

Automatic shift to lateral, AP and oblique angles (only Xmplar-dr)

Machine Dimensions 1500kg/3306lbs; LxWxH: 2810mm x 1586mm x 2276mm (111" x 63" x 90")

Diagnostic Viewing Station High Luminance, High Contrast, 23" Monitor; 1920 x 1080 pixels

Imaging software allowing imaging recall, database access, zoom, pan, rotate, lucid™ image

enhancement, anatomical measurement capability and window/level control.

Ambient Environment 10°C - 25°C, 40% - 75% relative humidity, 100hPa - 1060hPa atmospheric pressure, dust free

non-corrosive atmosphere.

Compliance US FDA Pre-Market Clearance (510k)

No. K013999 Accession # 0310920 CE-Mark, ISO 13485-2003, ISO 9001:2000 QMS Conformity Certificate # 5080212







C € 0120

info@lodox.com

7 Dartfield Road, Kramerville, Sandton 2146, South Africa. tel +27 11 444 9118 Lodox NA LLC, 143 Burton Street, Painesville, Ohio 44077. tel (440) 209-1246

www.lodox.com







MKT-19-0016-A/AUG2014

These specifications are current and subject to change. Lodox Systems reserves the right to revise these specifications to ensure a well engineered product.



