

# Measure Breath Nitric Oxide for Airway Inflammation with the NObreath<sup>®</sup> FeNO monitor



CE  
0086



*Improving asthma management, one breath at a time.*

[www.bedfont.com](http://www.bedfont.com)



# NObreath® FeNO monitor

## Fractional Exhaled Nitric Oxide (FeNO):

Airway inflammation is a central process in asthma and other lung diseases<sup>1</sup>. Being able to measure this inflammation and monitor the patient's reaction to medication is regarded as a gold standard in the management of respiratory diseases. It is increasingly recognised that the measurement of FeNO in particular constitutes a novel way to monitor separate aspects of diseases. These include asthma, COPD and interstitial lung diseases, that are not assessed by other means, such as lung function<sup>2</sup>.

FeNO measurements to evaluate airway inflammation in asthma represents a significant advance in respiratory medicine<sup>3</sup>, but until now this has been an expensive test to deliver in everyday practice.

## Benefits of performing FeNO tests:

- Non invasive, quick and easy to perform<sup>3</sup>
- Shows patient's response to treatment, enabling the correct prescription of medication
- Shows patient compliance
- Aids in identifying patients who do/do not require ongoing treatment<sup>4</sup>
- Shown to be superior to the majority of conventional tests of lung function, such as peak flow recording and spirometry<sup>3</sup>





## Each NObreath® comes complete with:

**50 Mouthpieces** – specifically designed with the latest bacterial filtration to remove 99.9% of airborne bacteria from the patient's breath<sup>5</sup>. These can be used up to 3 times per patient, dramatically reducing the cost of testing for FeNO.



| Order Code | Description                       |
|------------|-----------------------------------|
| NTK50      | 50 Mouthpieces and 1 NObreathFlo™ |

**NObreathFlo™** – eye level flow indicator makes keeping a constant flow during exhalation easy, even for young children. Using the NObreathFlo™ and mouthpieces provided allows the user to comply with ATS/ERS guidelines for FeNO testing<sup>2</sup>.

**Carry case** – for protecting the NObreath® whilst in storage or when being transported.

**25 alcohol free Cleaning Wipes** – Wipes/gels containing alcohol cannot be used on the NObreath® or any of its components. **Order code – WIPES**



## NObreath® Features

### Battery Indicator

The NObreath® requires 3 AA batteries, making it totally portable and easy to transport.

### Colour Touch Screen

For quick and easy use, with visual prompts for patients whilst taking a test to ensure correct results every time.

### Adult and child profiles

to ensure the best sample times.

### Internal pump and NO scrubber

Enabling warm up and recovery time to be a maximum of 60 seconds by constantly presenting the sensor with NO free ambient air.



Measuring FeNO with NO<sup>breath</sup>® is as easy as  
1, 2, 3

1



**Inhale**

2



**Exhale**

3



**Readings  
instantly  
available**

# Technical Specification

|                                              |                                                                                                                                  |
|----------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------|
| Concentration range:                         | 5-300ppb nitric oxide                                                                                                            |
| Accuracy:                                    | ± 5ppb of measured value ≤50ppb<br>± 10% of measured value >50ppb                                                                |
| Repeatability:                               | ± 5ppb of measured value ≤50ppb<br>± 10% of measured value >50ppb                                                                |
| Sensor sensitivity:                          | 5ppb                                                                                                                             |
| Breath test time:                            | Adult 12 seconds<br>Child 10 seconds                                                                                             |
| Response time                                | <10 secs                                                                                                                         |
| Warm up time                                 | ~60 secs                                                                                                                         |
| Ambient air test:                            | 30 seconds                                                                                                                       |
| Operating temperature range:                 | 10-30°C (ambient)                                                                                                                |
| Operating relative humidity (environmental): | 10-80% Rh (non-condensing)                                                                                                       |
| Sensor operating life:                       | 1-2 years; 6 month warranty                                                                                                      |
| Detection principle:                         | Electrochemical sensor                                                                                                           |
| Sensitivity drift:                           | <5% per annum                                                                                                                    |
| Maximum ambient operating level:             | 350 ppb NO                                                                                                                       |
| Power:                                       | 4.5V DC: 3 x AA (LR6 or equivalent) alkaline batteries                                                                           |
| Battery Life (1 set of 3 AA batteries):      | Up to 120 tests                                                                                                                  |
| Display:                                     | Colour LCD with touch screen                                                                                                     |
| Dimensions:                                  | Approx. 152 x 87 x 47mm                                                                                                          |
| Weight:                                      | Approx. 400g including batteries                                                                                                 |
| Construction:                                | Case – Polycarbonate/ABS blend with elastomeric overmould<br>NObreathFlo – Polycarbonate/ABS blend<br>Mouthpiece – Polypropylene |

<sup>1</sup>Shelhamer JH, Levine SJ, Wu T, Jacoby DB, Kaliner MA, Rennard SI. NIH conference: airway inflammation. Ann Intern Med 1995;123:288-304.

<sup>2</sup>ATS/ERS Recommendations for Standardized Procedures for the Online and Offline Measurement of Exhaled Lower Respiratory Nitric Oxide and Nasal Nitric Oxide, 2005; American Journal of Respiratory and Critical Care Medicine; vol. 171: 912-930;2005

<sup>3</sup>Andrew D. Smith, Jan O. Cowan, Sue Filsell, Chris MacLachlan, Gabrielle Monti-Sheehan, Pamela Jackson and D. Robin Taylor. Diagnosing Asthma: Comparisons between Exhaled Nitric Oxide Measurements and Conventional Tests. Am J Respir Crit Care Med Vol 169. pp 473-478, 2004.

<sup>4</sup>D R Taylor, MW Pinenburg, A D Smith and J C D Jongste. Exhaled nitric oxide measurements: clinical application and interpretation. Thorax 2006;61:817-827.

<sup>5</sup>Health Protection Agency (HPA). Porton Down, Report No. 43/06. Pp. 10-11



Contact Bedfont or one of our worldwide **NObreath®** distributors for a free demonstration

[www.bedfont.com](http://www.bedfont.com)  
Tel: +44 (0)1622 851122  
E-mail: [ask@bedfont.com](mailto:ask@bedfont.com)

A full list of our worldwide distributors can be found at  
[www.bedfont.com/distributors/](http://www.bedfont.com/distributors/)

***Scientific contributions to health.***

**Bedfont Scientific Ltd**  
Station Road, Harrietsham, Maidstone,  
Kent, ME17 1JA, England  
Tel: +44 (0)1622 851122, Fax: +44 (0)1622 854860  
Email: [ask@bedfont.com](mailto:ask@bedfont.com)  
[www.bedfont.com](http://www.bedfont.com)



© Bedfont Scientific Limited 2011

Issue 9 - April 2016, Part No: MKT140  
Bedfont Scientific Limited reserve the right to change or update this literature without prior notice.  
Registered in: England and Wales, Registered No: 1289798