




Always connected wireless,
with extremely low battery consumption.

Smart One[®]

Peak Flow and FEV1
now available on your Smart Phone

Simple and lightweight to carry on the go
Smart One[®] has been designed to satisfy the growing demand
for wireless devices capable to connect to the new generation
of Android and iOS Smart Phones and Tablets



MIR SMART ONE




MIR Smart One App includes
an **exclusive Incentive for adults
and children** based on both exhaled
Flow and Volume that is helpful
in improving subject's compliance



*Ideal for monitoring respiratory illnesses
and in the self management of Asthma,
COPD, Lung Transplant care, Cystic
Fibrosis and for use in Clinical Trials*



New turbine flowmeter
with plastic mouthpiece for personal use
easy to remove for cleaning



OEM Versions Available
MIR offers SDK API to allow third party
developers to run their own App



Graphic trend of Peak Flow or FEV1
Symptoms selection and easy scoring
system with the option of multiple
tapping to identify the intensity

*FDA clearance pending